



**Myo Munchee  
BEBE**  
6-18 mths

Suitable for all babies from 6 to 18 months old. If your baby is on the cusp of 18 months it is best to start with the Bebe and move on to the Mini after six months of Bebe use.

 **Please note**

that we only sell the Bebe Munchee if the patient is **under care of a Munchee trained practitioner.**



**Myo Munchee  
MINI**  
18 mths-4 yrs

Suitable for children from 18 months old and can include children up to 5 years old if their jaw width is narrow and they have a history of:  
-lip or tongue tie,  
-mouth breathing,  
-food avoidance, or  
-picky eating.



**Myo Munchee  
JUNIOR**  
4-8 yrs

Suitable for children from 4 to 8 years but can include 9 year olds who may have a mouth breathing habit, low oral tone, tooth crowding and/or lip/tongue tie which has impacted their jaw width and resulted in a more narrow jaw.



**Myo Munchee  
TWEEN**  
8-12 yrs & adults with  
narrow jaw structure

Designed for children from 8 years old with molar pads that can be removed for a comfortable fit if their 12 year old molars haven't come in.  
It may also be the most suitable size for an older teenager who has crowding and low oral tone, low tongue posture<sup>^</sup> and/or a history of lip/tongue tie.  
Transition from Tween to Adult is usually between 12-14 years.



**Myo Munchee  
ADULT**  
From 14 years+

Suitable for most adults and accommodates people from 12 to 14 years of age.  
Please refer to the Width Measurement Chart below, to confirm suitability for this size.  
**Size: Medium**  
**Size: Large\***  
  
\* for adults with broad skeletal structure

**Does your child have special needs?**

It is best to begin with a smaller size to make the introduction and management of the Munchee easier and more comfortable.

\*Low oral tone - Weak lips, tongue or cheeks, you may see this as dribbling/drooling.

<sup>^</sup>Low tongue posture - Is the tongue resting on the floor of the mouth? It should rest on the roof of the mouth.

**Myo Munchee Device Inter-Canine Width Measurements**

**BEBE**

N/a  
use age guide

**MINI**

N/a  
use age guide

**JUNIOR**

30 - 35mm

**TWEEN**

33 - 37mm

**ADULT**

Medium 35 - 40mm  
Large 40 - 45mm